
Lyset™

BOOST YOUR CELL CULTURE TODAY
FOR THE EXPERIMENTS OF TOMORROW

- platelet derived supplement of human origin
- fully standardized and lyophilized
- enables selection expansion and differentiation of human primary cell cultures
- allows culture of cells in a xeno-free environment

Lyset™ is manufactured
by Sclavo Diagnostics International S.r.l.
on Biorigen S.r.l. licence.

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Parmi les différentes alternatives au sérum animal, platelet derived preparations have been proposed since more than a decade. Nevertheless, the lack of standardized manufacturing and quality control procedures and the high degree of variability between the different platelet preparations, did not allow to reach a consensus on the applicability of this innovative cell culture medium supplement. Lyset™ is a platelet derivative obtained starting from human certified buffy coat samples with a defined platelet concentration and following standardized protocols including also lyophilization and biological activity testing, prior the product release as cell culture medium additive.

Lyset™ is being utilized to replace animal and human serum such as Fetal Bovine Serum (FBS) in culture medium formulation for different types of cells. The Lyset™ was investigated in culture studies (cell growth, viability and product release) towards a number of target cells including fetal and adult stem cells, articular chondrocytes, osteoblasts, skin fibroblasts and other primary cells. In general, Lyset™ supplemented medium supported cell growth and maintained viabilities superior to fetal bovine serum supplemented medium.

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Method

The cell supplement is obtained by the combination of two components: Lyset™ and Lyset Active Diluent™.

The lyophilized Lyset™ product restored with H_2O has a platelet equivalent concentration of about $1 \times 10^7/\mu$ l. In a standard preparation of the restored product the PDGF-BB and VEGF concentrations (quantified by ELISA assay) are at least 100 ng/ml for PDGF-BB and at least 2 ng/ml for VEGF. The lyophilized Lyset Active Diluent™ restored with H_2O has a platelet equivalent concentration of less than $5 \times 10^4/\mu$ l. The possibility of mixing the two supplement components in different percentages, allows the supplement customization and optimization for each cell type.

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
Lyset™ is particularly effective to increase the number of cells that can be derived from a tissue biopsy (primary cell cultures).

Cell doublings of human primary cells in culture. Osteoblasts; Adipose Derived Stem Cells; Umbilical Cord Derived Mesenchymal Stem Cells; Skin Fibroblasts.

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